


# Flexibility



Accurately goes where you need it.  
Stays where you place it.

**Artisan™**  
Control Catheter

**hansen®**  
MEDICAL  
The Global Leader in Flexible Robotics

# Flexibility Meets Stability.

## Artisan™ Control Catheter

The Artisan Control Catheter is a delivery system for percutaneous mapping catheters.\* Comprised of both inner and outer guides, the Artisan Catheter provides flexibility to reach difficult-to-access cardiac anatomy while maintaining stability during electrophysiology procedures.

The Artisan Catheter is driven by the Sensei™ Robotic Catheter System, which is designed to facilitate manipulation, positioning, and control of mapping catheters inside the heart.




## Artisan™ Control Catheter

\*Accommodates physician-preferred catheters 8F or smaller and a minimum of 115cm useable length.

### Outer Guide

### Inner Guide



Proximal Bend  
Max 20°

Distal Bend  
Max 90°

270° Deflection  
in Any Direction



800 East Middlefield Road | Mountain View, CA 94043 | PH: 650 404 5800 | FAX: 650 404 5901 | [www.hansenmedical.com](http://www.hansenmedical.com)

**The safety and effectiveness of the Sensei Robotic Catheter System for use with cardiac ablation catheters in the treatment of cardiac arrhythmias, including atrial fibrillation, have not been established.**

© 2008 Hansen Medical, Inc. All rights reserved. LT0001 091508  
"Artisan", "CoHesion", and "Sensei" are trademarks of Hansen Medical, Inc. and "Hansen Medical" and "Hansen Medical with heart logo" are registered trademarks of Hansen Medical, Inc., in the United States and other countries.